


TrichoScan training
Professional images
Camera: Fotofinder



TrichoScan

.....making hair growth measurable

Correct Image



Tattoo within the area

Black hairs (after hair dye)

No air bubbles

No dye remnants


No crossing hairs from outside

Sharp and bright image


Correct Image




Correct Image




Correct Image



Correct Image



Incorrect Image



Air bubbles

Dye remnants


Incorrect Image

Hair from the outside cross the
Measurement area




Incorrect Image

Dark spot in the middle





Incorrect Image

Air bubbles,
Image not in focus




Incorrect Image

Unclipped hairs from outside



Incorrect Image

Dye remnants




Incorrect Image

Image too dark, the ring light was
Switched off.

Incorrect Image




Hairs from outside




Incorrect Image

Air bubbles




Incorrect Image




No hair dye, wrong hair dye

Incorrect Image

Dye remants




Incorrect Image



Air bubbles


Incorrect Image

Dye remnants




Incorrect Image

Dye Remnants



Incorrect Image


Air bubbles



Incorrect Image


Unclipped hairs

Air bubbles



Incorrect Image


Dye Remnants



Incorrect Image

Unclipped hairs

Air bubbles





TrichoScan

.....making hair growth measurable